MULINO A TAGLIENTI SM 50

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MULINO A TAGLIENTI SM 50 IN CUSTOMER USE

Find out which processes our customers use the SM 50 for.

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The research group uses the SM 50 cutting mill to shred plastic waste, biomass, and electronic scrap for the development of new materials for a sustainable world.

"We use the SM 50 cutting mill to process plastic waste, biomass, and electronic scrap for the development of new materials that support a more sustainable future. The SM 50 allows us to quickly and easily pre-grind a wide variety of samples, enabling mechanochemical processes that ultimately lead to the creation of new materials in a circular mindset — a true asset to our scientific work."

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
MULINO A TAGLIENTI SM 50


FASTER AND MORE FLEXIBLE SAMPLE PREPARATION

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Different tasks, different users


In contract laboratories or academic institutions, many people use the same laboratory grinder for a wide variety of samples. It can happen that feathers are ground today and textiles or completely different sample types tomorrow.

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Optimized for different sample tasks in a row

Healthy nutrition is the basis for a happy life. Our food should combine enjoyment and health. To achieve this, it is important that our food is free from harmful substances and that we have precise knowledge of its composition. The quality of agricultural products is tested in environmental and food laboratories.

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Grinding of small sample quantities

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ARTICOLI

CUTTING MILL SM 50

Cutting Mill SM 50 (please order milling chamber set, rotor, hopper and bottom sieve separately)

20.720.0001  SM 50 200-230V 50/60Hz, with frequency converter for variable speed

20.720.0003  SM 50 110V 50/60Hz, with frequency converter for variable speed

MILLING CHAMBER SETS SM 50

Consisting of milling chamber, double acting cutting bars, door insert and 500 ml glass bottle / receptacle

22.010.0003 For standard applications: Aluminum milling chamber with stainless steel double acting cutting bars, stainless steel cyclone, stainless steel door insert, and 500 ml glass bottle

22.010.0004 For abrasion and corrosion resistant applications: Stainless steel milling chamber with stainless steel double acting cutting bars, stainless steel cyclone, stainless steel door insert, and 500 ml glass bottle

22.010.0005 For grinding without heavy metal contamination: Aluminum milling chamber with steel 1.1740 double acting cutting bars, stainless steel cyclone, heavy metal free steel door insert, and 500 ml glass bottle

22.010.0006 For grinding without cyclone: Aluminium gravity chamber with stainless steel double acting cutting bars, stainless steel door insert, adapter incl. stainless steel 3 l receptacle

ROTORS SM 50

22.608.0041 3-disc rotor made of stainless steel with reversible cutting tips made of tungsten carbide

22.608.0042 3-disc rotor made of steel 1.0570 with reversible cutting tips made of tungsten carbide, for grinding without heavy-metal contamination

22.608.0039 Parallel-section rotor made of stainless steel with cutting plates made of stainless steel

22.608.0040 Parallel-section rotor made of steel 1.0570 with cutting plates made of steel 1.1740, for grinding without heavy-metal contamination

HOPPER SM 50

22.785.0011

EasyInspect hopper made of stainless steel with PE plunger

BOTTOM SIEVES SM 50

BOTTOM SIEVES MADE OF STAINLESS STEEL 1.4404

03.647.0574

Fori trapezoidali, 0.25 mm

03.647.0575

Fori trapezoidali, 0.5 mm

03.647.0576

Fori trapezoidali, 0,75 mm

03.647.0577

Fori trapezoidali, 1,00 mm

03.647.0578

Fori quadrati, 2,00 mm

03.647.0579

Fori quadrati, 4,00 mm

03.647.0580

Fori quadrati, 6,00 mm

03.647.0581

Fori quadrati, 8,00 mm

03.647.0582

Fori quadrati, 10,00 mm

22.647.0001

Bottom sieves made of stainless steel complete set of 9 bottom sieves 0.25 - 10 mm

BOTTOM SIEVES MADE OF STEEL 1.0353 FOR GRINDING WITHOUT HEAVY-METAL CONTAMINATION

03.647.0583

Fori trapezoidali, 0.25 mm

03.647.0584

Fori trapezoidali, 0.5 mm

03.647.0585

Fori trapezoidali, 0,75 mm

03.647.0586

Fori trapezoidali, 1,00 mm

03.647.0587

Fori quadrati, 2,00 mm

03.647.0588

Fori quadrati, 4,00 mm

03.647.0589

Fori quadrati, 6,00 mm

03.647.0590	Fori quadrati, 8,00 mm
03.647.0591	Fori quadrati, 10,00 mm
22.647.0002	Bottom sieves made of steel 1.0353 complete set of 9 bottom sieves 0.25 - 10 mm

ACCESSORIES SM 50

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22.523.0001		Sample bottle 250 ml (10 pcs.)
22.523.0002		Sample bottle 500 ml (10 pcs.)
02.107.0493		Lid for collecting receptacle cyclone 3l (only for cyclone)
02.011.0031		Collecting receptacle 3 litres, stainless steel (for grinding with or without cyclone)
02.183.0146		Adapter for collecting receptacle 3l and 5l (only for SM 50 gravity outlet)
02.011.0029		Collecting receptacle 5 litres, stainless steel (only for SM 50 gravity outlet)

MILLING CHAMBERS

22.010.0007	Milling chamber for cyclone made of aluminium with stainless steel double acting cutting bars for standard grinding applications (A)
22.010.0008	Milling chamber for cyclone made of aluminium with steel 1.1740 double acting cutting bars for grinding without heavy metal contamination (A)
22.010.0009	Milling chamber for cyclone made of stainless steel with stainless steel double acting cutting bars for corrosion resistant applications (A)
22.010.0010	Milling chamber for gravity use made of aluminium with stainless steel double acting cutting bars (B)
22.730.0004	Door insert made of heavy metal free steel 1.0570
22.730.0003	Door insert made of stainless steel

STATIONARY CUTTING BARS FOR MILLING CHAMBERS

22.152.0012	Double acting cutting bars made of stainless steel (3 pcs.) for SM 50
-------------	---

22.152.0013	Double acting cutting bars made of steel 1.1740 (3 pcs.) for grinding without heavy metal contamination
-------------	---

22.608.0043	Adjustment tool for double acting cutting bars
-------------	--

ACCESSORIES FOR ROTORS

22.908.0001	 Reversible cutting tip incl. screws, tungsten carbide, for 3-disc rotor (10 pcs.)
-------------	---

22.151.0010	Cutting plates made of stainless steel for parallel section rotors (3 pcs.)
-------------	---

22.151.0011	Cutting plates made of steel 1.1740 for parallel section rotors (3 pcs.)
-------------	--

CYCLONES

22.935.0040	Cyclone made of stainless steel
-------------	---------------------------------

22.748.0005	Industrial vacuum cleaner for mills, 230 V, 50/ 60 Hz (other electrical versions and spare filter bags available upon request)
-------------	--

ADDITIONAL ITEMS

03.725.0047	PE plunger for EasyInspect hopper
-------------	-----------------------------------

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99.200.0045	IQ/OQ Documentation for SM 50
-------------	-------------------------------