LUCHTSTRAALZEEFMACHINE AS 200 JET

The Air Jet Sieving Machine AS 200 jet is particularly suitable for low density and low particle size materials which tend to agglomerate. It is used with sieves of 10 microns mesh size and up. The procedure is very gentle on the material as no mechanical sieving aids are required. The average sieving time is only 2-3 minutes.

The AS 200 jet is specifically designed for test sieves with a diameter of 203 mm/8" (or 200 mm with adapter). The air jet generated by an industrial vacuum cleaner can be adjusted by using the manual vacuum regulation. Optionally, an automatic vacuum regulation is available. The Open Mesh Function, a procedure which greatly reduces the number of near-mesh particles, provides optimum separation efficiency, excellent reproducibility and a longer service life of the sieves. Sieving time and nozzle speed are conveniently selected with a single button; the settings are shown in the graphic display. The Quick Start Mode is used to start the sieving process under standard conditions without entering parameters.

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- | Air jet technology for dispersion and deagglomeration
- | Measuring range 10 µm to 4 mm
- | Quick, efficient procedure
- | Open Mesh Function reduces near-mesh particles
- | Digital parameter setting (time, vacuum, speed)
- | Quick Start option
- | Variable nozzle speed
- | Automatic vacuum regulation and cyclone (options)
- | Memory for 9 Standard Operating Procedures (SOPs)
- | Suitable for RETSCH standard sieves

LUCHTSTRAALZEEFMACHINE AS 200 JET

TOEBEHOREN & OPTIES

RETSCH sieve shakers can be equipped with a variety of accessories to meet a wide range of application requirements.

Cyclone with holder and collecting receptable

To extend the service life of the filters in the vacuum cleaner and for recovery of the sample fraction passing the sieve, we recommend the use of the optional cyclone. The separation degree and limiting particle size respectively depend on the sample characteristics.

Test sieves 10 µm and up

RETSCH test sieves consist of a solid stainless steel sieve frame of high stability for reliable sieving results. Paying close attention to mesh-specific requirements, the fabric of each test sieve is precisely joined into the frame and tautened.

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RETSCH TEST SIEVES AND ACCESSORIES - ENGINEERED FOR SUPERIOR PERFORMANCE

LUCHTSTRAALZEEFMACHINE AS 200 JET

TYPISCHE MONSTER MATERIALEN

RETSCH's Air Jet Sieve AS 200 jet is ideally suited for separation, fractioning and particle size determination of pharmaceutical materials, ceramics, chemical products, cosmetics, food, minerals, pigments, plastics, powder coating, rubber, toner particles and washing powder.

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

EASYSIEVE / EASYSIEVE CFR EVALUATIE SOFTWARE

EasySieve, the software for particle size analyses, exceeds manual evaluation in many aspects. The software is able to automatically control the necessary measurement and weighing procedures – from the registration of the weight of the sieve up to the evaluation of the data. It is simple and convenient to use and is also available in an FDA 21 CFR Part 11-conform version.

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BESTELGEGEVENS

LUCHTSTRAALZEEFMACHINE AS 200 JET

Air Jet Sieving Machine for test sieves 203 mm (8") Ø, with manual vacuum regulation, incl. lids for test sieves 1" and 2" height and rubber mallet
(please order test sieve and vacuum cleaner separately)

30.027.0001		AS 200 jet	100–240 V, 50/60 Hz	incl. test rapport vlgs. EN 10204 2.2
-------------	---	------------	---------------------	---------------------------------------

VACUUM CLEANER AS 200 JET

32.748.0004		Industrial vacuum cleaner GM 80, 230 V, 50/60 Hz (other electrical versions available) only for 30.027.0001
-------------	---	---

32.524.0010		Main filter for industrial vacuum cleaner GM 80, 1 piece
-------------	---	--

32.524.0011		HEPA filter unit for industrial vacuum cleaner GM 80, incl. 1 HEPA filter
-------------	---	---

32.524.0012		HEPA filter for industrial vacuum cleaner GM 80, 1 piece
-------------	---	--

32.100.0004		Automatic vacuum regulation AS 200 jet
-------------	--	--

CYCLONE AS 200 JET

32.935.0008		Cyclone with holder and collecting receptacle (250 ml and 500 ml sample bottle) (industrial vacuum cleaner GM 80 (32.748.0004) recommended) only for 30.027.0001
-------------	--	---

ACCESSORIES AS 200 JET / JET PRO / JET PHARMA

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22.523.0002

Monsterflessen, 500 ml, 10 stuks

TEST SIEVES WITH ELECTROFORMED SHEET

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ADDITIONAL ITEMS AIR JET SIEVING MACHINES

99.200.0003

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RECALIBRATION AIR JET SIEVING MACHINES

The Retsch Air Jet Sieving Machines can be calibrated at Retsch/Germany or at selected subsidiaries/agents

SER-97.001.0001

Calibration of AS 200 jet at Retsch/Germany
incl. test report acc. to EN 10204 2.2